

EDU™ Series



EDU™ Series Ductless Demonstration Fume Hoods

- Safe for Education, Vocational Training and Industrial Workforce Development
- Meets or Exceeds OSHA, ANSI and other International Standards



60 watt¹ EDU-MOBILE, EDU-M-60



24–62 watt¹

The single EC blower motor assures lower cost of ownership in one of the world's most energy efficient ductless fume hoods.



"The World's Most Extensive Selection of Containment Solutions."



CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)

Get a Quote. 

INTRODUCTION

The EDU Series Ductless Demonstration Fume Hoods are designed to provide 360° visibility while protecting users and the classroom environment from hazardous vapors generated on and above the work surface.



24–62 watt¹

The single EC blower motor assures lower cost of ownership in one of the world's most energy efficient ductless fume hoods.

APPLICATIONS

Using innovative filtration technology, the EDU Series creates a safe work environment over the widest range of applications in the industry.

Science Class Lectures / Safety Practice
Demonstrations / Laboratory Assignments



Deep into its second generation, Air Science embraces the diversity and cultural heritage of the founders and co-workers who are continuing a tradition of excellence. Demonstrating a commitment to adaptation, inclusion, and quality output from a United States-based company with a domestic and global reach.



KEY FEATURES

- High efficiency ebm-papst EC blower.
- Energy saving LED lighting.
- Protects the operator from fume and (with added HEPA/ULPA filter) particle hazards.
- Improved filter clamping eliminates bypass leakage.
- Low airflow alarm.

PERFORMANCE ADVANTAGE

Advanced carbon filtration technology offers a safe, high performance alternative to conventional ducted fume hoods for a broad range of applications.

Environmental Benefits. Air Science® ductless fume hoods isolate and trap chemical vapors to prevent ecological impact through release into the environment.

Versatile. Each filtration system is selected for its specific application. Carbon filters are available in more than 14 configurations for use with vapors of organic solvents, acids, mercury, and formaldehyde. HEPA/ULPA filters can be added for biological safety.

Easy to Install. The ductless fume hood is self-contained and does not require venting to the outside. Many units are portable and may be moved with minimal downtime and without filter changes. Set-up, operation, and filter maintenance are straightforward.

Energy Efficient. Because filtered air is returned to the room, no demands are required of the facility HVAC capacity for make-up air.

Cost Effective. Facility ductwork, HVAC, and construction costs are eliminated.



39 watt¹ EDU-MOBILE, EDU-M-40

Safe to Use. Cabinet airflow and face velocity protect users from incidental exposures to fumes.

Self-Testing. (select models) Electronic airflow monitoring assures continuous safety. An electronic gas sensor monitors carbon filter performance.

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)

Get a Quote. 



39 watt¹ EDU-Mobile EDU-M-40, shown with optional duplex electrical outlet, gas petcock, cup sink and gooseneck water faucet.

DESIGN FEATURES

- A. Filter I.D. Window:** A convenient, strategically placed front cover window shows the installed filter part number and installation date to encourage timely filter replacement.
- B. Control Panel:** Electronic controls and displays include switches for the blower and low airflow alarm.
- C. Airflow Alarm:** Low airflow alarm continuously monitors filter loading and alerts user when service is needed.
- D. Air Velometer:** An analog air velocity meter is positioned in the user's field of vision.
- E. Steel Support Frame:** The chemical resistant epoxy coated steel frame adds mechanical strength. Optional all polypropylene construction is available if desired; see accessories.
- F. Hinged Front Sash:** When closed, the cabinet sash protects the contents from inadvertent external contact and better isolates the air within. The sash is easy to open and close.
- G. Work Surface:** The internal work surface is fitted with a standard polypropylene spill tray (available in white and black). An optional stainless steel tray is also available, see accessories.
- H. Pass Through Ports:** Electrical cords and cables are safely routed into the cabinet through ports on the back.
- I. Electrostatic Pre-Filter:** The electrostatic pre-filter is accessible from inside the chamber and 91% effective down to 1-3 microns.
- J. Filter Door Key:** Filter access keys prevent unauthorized removal or accidental exposure to dirty filters.

- K. Internal Manual Speed Controller:** Authorized personnel set the centrifugal fan motor speed as desired.
- L. Locking Casters:** Permit complete portability from lab to lab. The overall width and height are sized to fit easily through standard doorways.
- M. Isolated Single Switch Electrical Sockets:** 13 amp.
- N. Quick release connections:** For water, waste, and gas.
- O. Optional petcocks are available for gas and water.** An optional gooseneck faucet is positioned over an optional chemically-resistant polypropylene sink cup integrated into the work surface.

ADDITIONAL FEATURES

360 Degree Visibility: Clear back and side panels allow ambient illumination into the chamber and provide users with an unobstructed view of its contents.

Construction: All models are available in either metal or polypropylene construction. See selection chart for specifications and dimensions. Available in 120V, 60Hz and 230V, 50Hz models.

No Installation, No Ducting Required: Self-contained, integrated systems ideal for fixed location or mobile applications.

Safe, Energy-Efficient Operation: All conditioned air is safely returned to the room, minimizing impact on facility HVAC costs.

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)

Get a Quote. 

PERFORMANCE

Air Science [Multiplex Filtration](#) offers a range of options for high performance protection.

Multiplex filter configuration permits a customized combination of filter media for a broad range of chemical families and biological agents if required.

EFT™ filtration technology broadens the Air Science application for ductless fume hoods.

A high capacity air handling system delivers face velocity of 100 fpm.

DESIGN

Professional quality Air Science fume hoods comply with current technical and safety regulations. The cabinet frame and work surfaces, comprised of industrial components, are durable and chemically resistant.

The Air Science filter assembly is easy to access, easy to change, plus a unique filter clamping design eliminates bypass leakage outside the cabinet.

The optional SafeSwitch HEPA Filter Shutter System is available for safer filter exchange.

RELIABILITY

Internal systems are isolated from fumes, extending product life.











Energy-efficient ebm-papst brand EC blowers promote long life and dependable performance of EDU Series fume hoods.

SELECTION

EDU Series products are available in 4 configurations and 6 standard sizes, in metal or polypropylene construction, totaling 12 standard models.

EDU SERIES DIFFERENTIAL SUMMARY

				
Feature	EDU-MOBILE	EDU-CLASSIC	EDU-ADA	EDU-JUNIOR
General Description	Mounted on wheeled cart with small side storage cabinet. Blowers, filters and controls are mounted below the work surface to lower the center of gravity.	Mounted on a wheeled cart with an enclosed chemical storage area in the base cabinet.	Provides wheelchair access to the workstation front. Controls are mounted in a handheld remote box placed inside the workzone.	A basic workstation, mounted on wheeled cart.
Airflow				
Air Sampling Port	•	•	•	
Base Cabinet	•	•		
Dynamic Filter Chamber	•	•	•	
Airflow Alarm	•	•	•	•
Hour Counter	•	•	•	
Optional Service Fixtures, Cup Sink	•	•	•	•
Push/Pull Handle	•	•		
Dwyer Continuous Airflow Display	•	•	•	•
Safety Filter	•	•	•	
Storage Compartment	•	•		
Track and Wheel Filter Insert/Remove Function	•	•	•	

CONTENTS:


- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)

Get a Quote. 



DOCKING STATION

An optional docking station prevents accidental disconnect of services to mobile fume hoods and prevents unauthorized access to main service connections which are securely locked away when not in use.



Basic Control Panel



Advanced Control Panel



FSA/Autocal Control Panel



Monitair Control Panel

CONTROL

The **basic control panel** is standard on the EDU-JUNIOR and includes an On/Off switch and low airflow alarm.

The **advanced control panel** includes an On/Off switch, low airflow alarm and hour meter to aid in determining available filter life. Standard on EDU-MOBILE, EDU-CLASSIC, and EDU-ADA.

The optional **FSA/Autocal controller** displays the airflow and uses an electronic gas sensor to detect when the filter needs changed. Audio and visual alarms alert users to filter saturation and if the airflow reaches preset thresholds.

The optional **Monitair microprocessor controller** monitors and displays cabinet operating parameters, airflow, containment, and filter condition; emits audio and visual alerts if conditions become unsafe, all on an LCD display.

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)

Get a Quote. 



FILTRATION

At the heart of the EDU Series product line is innovative filtration technology. **The Multiplex Filtration System** consists of a pre-filter, main activated carbon or HEPA/ULPA filter, and safety activated carbon or HEPA/ULPA filter. The system permits a customized combination of filter media and configuration for chemical and physical adsorption specific to each application need.

The Air Science **carbon filtration technique** is based on enhanced, activated carbon particle formulations from specially selected, naturally occurring raw material that is superior to wood or other organic sources. The carbon is treated to attain the proper porosity and aggregate surface area and to react with several ranges of aerosolized chemicals moved through the filter by an air handling blower.

View available filters and descriptions on [page 11](#).



safe SWITCH™
The optional SafeSwitch HEPA Filter Shutter system ensures that operators are safely separated from trapped contaminants during filter changes.

FILTER CONFIGURATION

The Multiplex feature permits one or more filtration options to be combined to meet a wider range of multiple-use applications.

The EDU Fume Hood can be equipped with a single activated carbon main filter or with a stacked configuration which combines two main filters, each activated to adsorb one or more specific vapors or family of vapors. For safety against particulates, an optional HEPA or ULPA can also be added. When used with a HEPA/ULPA filter, the ductless fume hood may be applied as a Class I Biological Safety Cabinet.

The carbon filter is sized to fit the specified product model number and configured to optimize airflow across 100% of the filter surface area. The self-contained assembly maximizes filter efficiency, prolongs filter life, optimizes diffusion and saturation, and improves user safety.

- P. Electrostatic Pre-Filter:** Protects the main filters from aerosols, mists, dust, and particulates.
- C. Activated Carbon Main Filter:** A single, blended, or stacked filter configuration.
- H. HEPA/ULPA Filter, Optional:** Both HEPA and ULPA filters use micro-glass fiber media designed to capture fine particles and biologicals. Both filters can capture particles smaller than the micron size for which they are tested. HEPA and ULPA filter efficiencies are 99.995% at 0.3 microns and 99.9995% at 0.12 microns respectively.

MULTIPLEX FILTRATION SYSTEM, SUMMARY


Application	Chemical	Powder/ Biological	Chemical & Powder	Chemical within Cleanroom
Safety Filter, Optional	C	H	H or C	H
Primary Filter	C	H	H C	H C
Pre-Filter	P	P	P	P

The system can be configured for the capture of acids, bases, and particulates, such as biological aerosols, when paired with HEPA or ULPA filters.




62 watt¹
EDU-ADA, P10XL-CART-RC-P

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
-  Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)

Get a Quote. 




EDU-MOBILE
DOWNFLOW

AIRFLOW

EDU Series fume hoods maintain a constant face velocity of 100 FPM in compliance with U.S.A. and International standards for safety and performance. Contaminated air is pulled through the Multiplex filtration system where activated carbon adsorbs chemical vapors, returning clean air to the room.

The main filters are easy to replace; no tools required. The filter clamps tightly against the filter gasket to prevent filter bypass and to maintain filter integrity.



EDU-CLASSIC, EDU-ADA and EDU-JUNIOR
UPFLOW

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)

Get a Quote. 




ENHANCED FILTRATION

The Air Science Enhanced Filtration Technology (EFT) is a universal filtration system developed for use with a wide range of core chemical families. These include organic acids, alcohols, aliphatic hydrocarbons, aromatic hydrocarbons, esters, aldehydes, ketones, ethers, halogens and others. Although the EFT system is weighted to accommodate these families, it can handle inorganic acids as well.

The Air Science EFT system is available as an option on Air Science EDU series ductless fume hoods.

Independent Test Results Independent testing confirms that the Air Science EFT system is superior in critical areas to other “green” fume hood systems recently introduced to the industry. AFNOR NFX 15-211 requires that three chemicals (isopropanol, cyclohexane, and hydrochloric acid) be tested under very precise conditions to ascertain and establish retention capacity at 1% of the threshold limit value (TLV) for each chemical.

Retention capacity (grams) for a single module at 1% of the TLV (Threshold Limit Value)

Specification	AFNOR NFX 15-211	
	IBR	Intertek
Testing Laboratory		
Product Manufacturer	Air Science	Brand E
Filter Type		Green
Test Results		
Isopropanol (alcohol)	2052	673
Cyclohexane (aliphatic hydrocarbon)	1531	914
Hydrochloric acid (inorganic acid)*	1205	2729*

**Based on “core” chemical families typically used in ductless fume hood applications, the Air Science EFT filter offers significant advantages over filters marketed as “universal” filters. With moderate to heavy acid applications, all ductless fume hoods made of metal are subject to corrosion and rust. On inorganic acids, the EFT filter provides a lesser, but more realistic, usable capacity.*



Avoid Revolving Filters Air Science strongly discourages the unsafe practice of revolving secondary backup filters into the primary filter compartment. All Air Science units are designed to avoid this false sense of security.

In a revolving filter system, users are instructed to rotate the secondary backup filter into the primary filter position after non-permissible exposure levels of chemicals are detected within the monitoring chamber.

Depending on when the unit can be properly shut down, the secondary filter can be loaded to the point of saturation itself, thereby creating a safety hazard if the filter is considered new.

If a new spare filter is not immediately available, a user may inadvertently (or knowingly) re-install a contaminated primary filter into the secondary location permitting the system to operate without protection.

Additionally, the secondary filter can become contaminated as it ages, sometimes for years, on top of an operational cabinet, losing filter efficiency by the time it is installed.

Either practice puts both personnel and the environment at risk, even though some manufacturers provide stickers to label the filters as “used.”

The Air Science non-revolving filter practice ensures that only a new filter is fitted into the primary filter compartment and permits the secondary filter to remain installed for at least twice the change-out period, resulting in a 50% savings in filter change-out costs.



secur.
safe disposal service



Filter disposal services are available in selected markets providing responsible destruction or recycling of saturated filters in authorized facilities.

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)

Get a Quote.




EDU-MOBILE


The EDU-MOBILE is mounted on a wheeled cart with a small side storage compartment. Blowers, filters and controls are mounted beneath the work surface to lower the center of gravity and simplify moving from room to room.

EDU-MOBILE, EDU-M-40


EDU-M-40




EDU-M-48



EDU-M-60



SIDEVIEW



MODEL	DIMENSIONS			WEIGHT (LBS/KG)	
	Internal Height	External (W × D × H)	Shipping (W × D × H)	Net	Ship
EDU-MOBILE					
EDU-M-40	38" / 965 mm	39.5" × 27.375" × 78" / 1003 × 695 × 1981 mm	48" × 40" × 90" / 1219 × 1016 × 2286 mm	270 / 123	320 / 145
EDU-M-48	38" / 965 mm	49.5" × 27.375" × 78" / 1250 × 695 × 1981 mm	60" × 40" × 90" / 1524 × 1016 × 2286 mm	370 / 168	420 / 191
EDU-M-60	38" / 965 mm	59.5" × 27.375" × 78" / 1511 × 695 × 1981 mm	72" × 40" × 90" / 1829 × 1016 × 2286 mm	470 / 213	520 / 236

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)


Get a Quote.




EDU-CLASSIC

The EDU-CLASSIC has horizontal inflow with top mounted filters and exhaust. This model is mounted on a wheeled cart with an enclosed chemical storage area in the cabinet base.

P15-XT-CART-MCC-40



SIDE VIEW



EDU-CLASSIC, P15-XT-CART-MCC-40

MODEL	DIMENSIONS			WEIGHT (LBS/KG)	
	Internal Height	External (W × D × H)	Shipping (W × D × H)	Net	Ship
EDU-CLASSIC					
P15-XT-CART-MCC-40	38" / 965 mm	39.5" × 28.5" × 87" / 1003 × 724 × 2210 mm	48" × 40" × 93" / 1219 × 1016 × 2362 mm	350 / 159	425 / 193

CONTENTS:


- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)

Get a Quote.



EDU-ADA

The EDU-ADA permits wheelchair access to the workstation. Controls are mounted in a handheld remote box that can be placed adjacent to or inside the work area.



EDU-ADA, P10XL-CART-RC-P

MODEL	DIMENSIONS			WEIGHT (LBS/KG)	
	Internal Height	External (W × D × H)	Shipping (W × D × H)	Net	Ship
EDU-ADA					
P10XL-CART-RC-P	38" / 965 mm	34" × 27.375" × 87" / 864 × 695 × 2210 mm	48" × 40" × 90" / 1219 × 1016 × 2286 mm	234 / 107	284 / 129

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)


Get a Quote.



EDU-JUNIOR

The EDU-JUNIOR is a basic workstation mounted on a wheeled cart.

P5-24XT-CART SIDE VIEW



EDU-JUNIOR, P5-24XT-CART

MODEL	DIMENSIONS			WEIGHT (LBS/KG)	
	Internal Height	External (W × D × H)	Shipping (W × D × H)	Net	Ship
EDU-JUNIOR					
P5-24XT-CART	24" / 610 mm	24" × 27" × 70" / 610 × 686 × 1778 mm	40" × 40" × 80" / 1016 × 1016 × 2032 mm	176 / 80	225 / 102

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)

Get a Quote. 

Warranty Info.

PRODUCT SPECIFICATIONS

EDU Series	EDU-MOBILE			EDU-CLASSIC	EDU-ADA	EDU-JUNIOR
Filtration	EDU-M-40	EDU-M-48	EDU-M-60	P15-XT-CART-MCC-40	P10XL-CART-RC-P	P5-24XT-CART
Airflow ¹	220 cfm	220 cfm	220 cfm	220 cfm	145 cfm	135 cfm
Pre-Filter	100 fpm	100 fpm	100 fpm	100 fpm	100 fpm	100 fpm

Construction

Construction	<... White epoxy steel frame and head unit, Clear sides and back panel. ...>					
Blower	<... ebm-papst EC blower. ...>					
Controls	<... Main On/Off. ...>					
Electrical	<... 120V, 60Hz or 230V, 50Hz voltages available. Specify when ordering. Other voltage options available. ...>					
Monitoring	<... Low airflow alarm, standard. ...>					

Efficiency

Power Consumption ¹	39 watt	59 watt	60 watt	39 watt	62 watt	24 watt
Lighting	<... LED lighting. ...>					
Noise, dBA ³	< 54	< 58	< 58	< 55	< 57	< 54

¹ All measurements are with Filter Type ASTM-030.

² Average airflow measured 6" (150 mm) from filter face.

FILTER SPECIFICATIONS

Purair Model	EDU-M-40	EDU-M-48	EDU-M-60	P15-XT-CART-MCC-40	P10XL-CART-RC-P	P5-24XT-CART
Safety Filter, Optional*	(1)	(1)	(1)	(1)	(1)	(1)
Primary Filter(s)*	(1)	(2)	(2)	(1)	(1)	(1)
Pre-Filter*	(1)	(1)	(1)	(1)	(1)	(1)

* For specific examples refer to Multiplex filtration system summary on [page 6](#).

FILTER SUMMARY

Formula	Description
GP Plus!	The most widely used filter in the range, primarily for solvent, organic, and alcohol removal.
ACI Plus!	Neutralizes volatile inorganic acid vapors.
ACR	Iodine and methyl iodide vapors as well as low level radioactive iodine.
ACM	Mercury vapor.
AMM	Removes vapors from dilute ammonia solutions and to remove low molecular weight amines.
SUL	Designed to remove hydrogen sulphide and low molecular weight mercaptans.
CYN	Removal of hydrogen cyanide. Many cyanide compounds will evolve HCN gas if acidified, so this filter is normally specified if working with any cyanide compound.
FOR	Designed to oxidize formaldehyde and glutaraldehyde fumes. It is widely used in hospital pathology laboratories.
EDU	Designed to handle chemicals normally used in a university level chemistry curriculum.
MIL	Designed for military applications involving war gasses.
HEPA/UPLA	Powders, particulates, and biologicals.
	Universal filtration.


CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)

Get a Quote. 

Warranty Info.

OPTIONS & ACCESSORIES

		EDU-MOBILE			EDU-CLASSIC	EDU-ADA	EDU-JUNIOR
Air Science Model		EDU-M-40	EDU-M-48	EDU-M-60	P15-XT-CART-MCC-40	P10XL-XT-CART-RC-P	P5-24XT-CART
Safety Filter	An additional carbon, HEPA or ULPA safety filter exceeding ANSI/AIHA Z9.5 requirements can be installed after the main filter.	<... Safety filters for vapor or particulate protection are available for all models. ...> Contact Air Science for ordering information.					
Monitair Controller*	Microprocessor controller monitors cabinet operating parameters, airflow, containment, and filter condition; emits audio and visual alerts if conditions become unsafe. Not TUV compliant.	MON-P	MON-P	MON-P	MON-P	MON-P	MON-P
Filter Saturation Alarm*	An electronic gas sensor emits audio and visual alerts when the main filter needs to be changed.	FSA	FSA	FSA	FSA	FSA	FSA
Spill Tray (Stainless)	Slides out for easy cleaning	SS-TRAY	SS-TRAY	SS-TRAY	SS-TRAY	SS-TRAY	SS-TRAY
SafeSwitch HEPA Filter Shutter System	Minimizes exposure to filter contaminants when removing used HEPA filters for insertion of new filters.	ASTM-030-SS	ASTM-030-SS	ASTM-030-SS	ASTM-030-SS	ASTM-030-SS	ASTM-030-SS
Polypropylene Construction*	Cabinets are available in all polypropylene construction. Contact Air Science for information.	EDU-M-40-PP	EDU-M-48-PP	EDU-M-60-PP	P15-XT-CART-MCC-40-PP	P10XL-XT-CART-RC-P-PP	P5-24XT-CART-PP
Duplex Electrical Outlet*	Two NEMA 14-20R receptacles with ground fault interrupter. 110V service standard; international fixtures available.	AS-GFI	AS-GFI	AS-GFI	AS-GFI	AS-GFI	AS-GFI
Service Fitting*	Cabinets can be fitted with service fixtures such as petcocks, faucets or valves	<... SF-X. Specify service fitting type (faucet, valve, petcock) and location when ordering ...>					
Cup Sink, Mounts into Tray*	Polyethylene cup sink (3" x 5" x 9") is fitted into the base tray. Other sizes and materials available. Contact Air Science for ordering information.	SINK	SINK	SINK	SINK	SINK	SINK
Docking Station* 	Security cable is shorter than utility and power line. Prevents damage from accidental disconnect of utility services to mobile fume hoods. Can be recessed or flush mounted. Includes key lock, quick release connections for water, waste and gas. An optional RJ45 outlet is available. Specify outlet choice when ordering.	DOCK-S	DOCK-S	DOCK-S	DOCK-S	DOCK-S	DOCK-S
Sliding Sash*	Vertical sliding sash permits easier access to the work surface when the mobile fume hood is used in limited space areas. Specify when ordering.	EDU-M-SS	EDU-M-SS	EDU-M-SS	P15-XT-CART-MCC-40-SS	P10XL-XT-CART-RC-P-SS	P5-24XT-CART-SS
Dual Access*	Units can be fitted with dual access front and rear panels for applications where student access is required.	NA	NA	NA	P15-XT-CART-MCC-40-DA	P10XL-XT-CART-RC-P-DA	P5-24XT-CART-DA

* Factory installed; specify when ordering

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Airflow Technology (p.6)
- Specifications (p.9)
- Options & Accessories (p.12)

Get a Quote. 

Warranty Info.

WARRANTY



This product is protected by the Air Science **Legacy Lifetime Warranty™** which starts on the date of shipment from our factory. This limited warranty is the result of thousands of successful Air Science production applications in pharmaceutical, laboratory, forensic, industrial, and educational applications.

This warranty covers defects in materials and workmanship. Our liability under the Legacy Lifetime Warranty is, at our option, to repair or replace any defective parts of this equipment if you document that these parts were defective at the time we sold the product to you. Normal conditions apply.

For details visit the [Service section](#) of our website at www.airscience.com.

STANDARDS & COMPLIANCE

Quality Management Systems	ISO 9001:2015
Chemical Fume Containment	ANSI/ASHRAE 110 1995
Carbon Filter Efficiency	BS 7989-2001 AFNOR NFX 15-211
Biological Safety Filter Efficiency HEPA and ULPA	IEST-RP-CC007.1 IEST-RP-CC001-4 EN 1822
Electrical Safety	TUV (except for units equipped with Monitair controllers)
Product Design	ANSI Z 9.5-2003 ANSI Z 9.7-1998
OSHA, Occupational Safety and Health Administration	OSHA Standard -29 CRF, Safety and Health Regulations for General Industry, 1910.1450: Occupational exposure to hazardous chemicals in laboratories. Part B, definition, laboratory type hood. All Air Science products meet this definition.
Environment	ISO 14001:2015 ENERGY STAR® Partner
Education (UK)	CLEAPPS Instruction Approved (EDU) based on ASTM-100 filter.



120 6th Street \ Fort Myers, FL 33907
T. 239-489-0024 \ Toll Free. 800-306-0656 \ F. 800-306-0677
www.airscience.com

The information contained in this manual and the accompanying product are copyrighted and all rights are reserved by Air Science. Air Science reserves the right to make periodic minor design changes without obligation to notify any person or entity of such change.