Pur Flow





- Meets or Exceeds OSHA, ANSI and Other International Standards

Purair FLOW-36







Air Science



Purair® FLOW Series laminar flow cabinets are designed to protect the interior work zone from particulate contamination.

The Purair FLOW Series product line employs the Air Science exclusive Multiplex™ HEPA Filtration technology to sustain a contamination-free environment.

HEPA filtered air passes uniformly through the cabinet interior via vertical laminar flow. The airflow is oriented to prevent introduction of airborne particulates during normal use.

APPLICATIONS

Vertical laminar flow cabinets are intended for use in non-hazardous applications where user protection from biologicals or biohazardous byproducts is not required.

Mycology and Food Microbiology \ Plant and Mammalian Cell Culture \ Clinical Pharmacies and Hospitals \ Cleanrooms \ Semiconductor Assembly \ Pharmaceutical Production \ Aerospace Contamination Control \ Medical Device Assembly \ General Research Protocols







KEY FEATURES

- Provides sterile work zone for aseptic techniques.
- Air cleanliness meets and exceeds ISO Class 5.
- Energy saving LED lighting.
- · Vertical laminar flow with HEPA filtration.
- Large front opening provides unrestricted access to the work zone.

FLOW SERIES

The Purair FLOW is designed to provide protection to the work surface from biological or particulate contamination. It is ideal for pharmacy compounding of non-toxic agents, assembly of finely tuned instrumentation, cell culture and other processes where protection of the user is not required but easy access to the work surface is important. ISO Class 5 air quality is maintained through a smooth, downward laminar flow that provides gentle yet uniform airflow throughout the cabinet. The clean, simple, low maintenance design offers a comfortable workstation for a range of applications and may be bench mounted or positioned on a mobile bench.



Purair FLOW-48 shown with optional spill tray and mobile cart.



Deep into its second generation, Air Science embraces the diversity and cultural heritage of the founders and co-workers who are continuing a tradition of excellence. Demonstrating a commitment to adaptation, inclusion and quality output from a United States-based company with a domestic and global reach.



DESIGN FEATURES

- A. Control Panel: The control panel includes an On/Off
- **B.** Main Filter: HEPA filter with 99.99% efficiency
- **C.** Pre-Filter: Disposable polyester fiber pre-filter
- **D.** Blower Motor: ebm-papst centrifugal fan.
- **E.** Stand: Optional base stand converts to mobile cart
- F. Lighting: Compact LED cabinet lamp located away
- **G.** Pass Through Ports: Convenient rear-wall pass through
- H. Filter Door Lock: Prevents unauthorized removal or
- **I.** Filter Clamp: Wide knob filter clamps are conveniently
- J. Work Surface: Optional polypropylene spill tray,

ADDITIONAL FEATURES

360 Degree Visibility: Clear back and side panels allow ambient illumination into the chamber and provide users

Standards Compliant: Performance specifications and

Construction: All models are available in either metal

Air Science FLOW cabinets are expertly designed to meet specific applications and certified for quality construction. Standard features, options and accessories are developed purposefully to enhance user-friendliness.

PERFORMANCE

The Air Science Multiplex filter offers a range of options for high performance protection.

Multiplex filter configuration permits a customized combination of filter media for a broad range of chemical families and biological agents if required.

Purair FLOW cabinets maintain up to 0.45 m/s (90 fpm) airflow velocity, measured 6" (150 mm) from the filter face.

The HEPA filters are easy to replace; no tools required.


DESIGN

Professional quality Air Science cabinets comply with current technical and safety regulations.

The cabinet frame and work surfaces, comprised of industrial components, are durable and chemically resistant with a microbial powder finish coating.

The Air Science filter assembly is easy to access, easy to change, plus a unique filter clamping design eliminates bypass leakage of the main filter.

Stackable sections make cabinets highly portable and enable quick setup.



SELECTION

Purair FLOW cabinets are available in three model sizes with various options. Designed for desktop use, or for installation on an optional base stand or mobile cart.

RFIJABILITY

Purair FLOW cabinets incorporate energy-efficient ebm-papst centrifugal blowers for maximum operational savings, low noise and minimal vibration.

Air Science uses long-life HEPA filters without aluminum separators to increase filter efficiency, minimize the potential for leakage and increase filter life.

CONTROL

The basic control panel features an On/Off switch for simple operation of the blower and light.



Basic Control Panel



FLOW

Purair FLOW cabinets incorporate energy-efficient ebm-papst™ permanently lubricated direct drive centrifugal blowers for maximum operational savings, low noise and minimal vibration.



Purair FLOW-36, shown with optional mobile cart.

Performance & Selection (p.4)

Containment & Filtration (p.5)

Options & Accessories (p.8





FILTRATION

At the heart of the Purair product line is innovative filtration technology. **The Multiplex Filtration System** consists of a pre-filter and main filter. The mechanical design enhances safety, convenience and overall value.

- The disposable pre-filter is accessible from the exterior top of the cabinet.
- A filter clamping mechanism allows for the filter to be easily installed and ensures an even seal at the filter peripheral face at all times to prevent bypass leakage.
- The filter chamber prevents contaminated air from contacting internal cabinet mechanisms.







Filter disposal services are available in selected markets providing responsible destruction or recycling of used saturated filters in authorized facilities.

AIRFLOW

Room air enters from the top of the cabinet through the disposable pre-filter where larger particles are trapped, increasing the service life of the main HEPA filter.

Air is forced evenly across the HEPA filter to deliver a flow of pure, uniform air within the work zone to dilute and flush airborne contaminants from the work area. A nominal filter face velocity of up to 0.45 m/s (90 fpm) ensures that there is a sufficient number of air changes within the work zone to eliminate cross contamination and to maintain optimum cleanliness.


Purified air travels across the work zone to the work surface in a vertical, unidirectional downflow stream and then exits the work zone across the open cabinet front.

FILTER CONFIGURATION

P. A The pre-filter may be replaced while unit is in operation.

H. The main filter is easy to replace; no tools required. The filter clamps tightly against the filter gasket to prevent filter bypass and to maintain filter integrity.

MULTIPLEX FILTRATION SYSTEM, SUMMARY		
Pre-Filter	P	Disposable polyester fibers with 85% arrestance.
Main Filter	H	A self-contained filter designed to physically capture particles larger than 0.3 microns (HEPA) or 0.12 microns (ULPA).






Through our partner company <u>Filtco Filters</u>, Air Science is a single source supplier of all pre-filters, carbon filters and HEPA/ULPA filters used in our products and those of many other manufacturers.


6

Specifications (p.6)
Options & Accessories (p.8)


Purair FLOW-24








Purair FLOW-48







MODEL	DIMENSIONS			WEIGHT (LBS/KG)			
	Nominal Width	Internal Height	Internal Depth	External (W x D x H)	Shipping (W \times D \times H)	Net	Ship
FLOW-24	24" / 610 mm	24" / 610 mm	24" / 610 mm	24" x 24" x 35" / 610 x 610 x 889 mm	40" x 40" x 45" / 1016 x 1016 x 1143 mm	72 / 33	150 / 68
FLOW-36	36" / 914 mm	24" / 610 mm	24" / 610 mm	36" x 24" x 35" / 914 x 610 x 889 mm	48" x 40" x 45" / 1219 x 1016 x 1143 mm	99 / 45	175 / 79
FLOW-48	48" / 1219 mm	24" / 610 mm	24" / 610 mm	48" x 24" x 35" / 1219 x 610 x 889 mm	60" x 40" x 45" / 1524 x 1016 x 1143 mm	138 / 63	225 / 102

Specifications (p.6)
Ontions & Accessories (p.8)

PRODUCT SPECIFICATIONS

Filtration	FLOW-24	FLOW-36	FLOW-48	
Airflow	< Vertical downflow; 0.3 m/s - 0.45 m/s (60 - 90 fpm)>			
Pre-Filter	< Disposable polyester fibers with 85% arrestance>			
Main Filter	<··· HEPA efficiency, 99.99% at 0.3μm. ···>			
Clamping	<··· Screw compression clamp. ···>			
Construction	FLOW-24	FLOW-36	FLOW-48	
Finish (exterior)	$<\cdots$ White epoxy-coated steel frame with white legs on cabinet sides. $\cdots>$			
Windows	< Acrylic, transparent>			
Blower	<··· ebm-papst centrifugal fan. ···>			
Controls	< Main On/Off switch for fan and lighting. Solid-state fan speed control with RFI filter maintains blower uniformity>			
Electrical	< 120V, 60Hz or 230V, 50Hz voltages available. Specify when ordering. Other voltage options available>			
Pass Through Ports	< Two, standard, knock-out>			
Efficiency	FLOW-24	FLOW-36	FLOW-48	
Lighting	<··· LED lighting. ···>			
Noise, dBA¹	< 63	< 66	< 66	

 $^{^{1)}}$ Measured 12" (30 cm) from the cabinet front and 15" (38 cm) above the work surface plane.

OPTIONS AND ACCESSORIES

Purair Model		FLOW-24	FLOW-36	FLOW-48
ULPA Filter	ULPA filter efficiency 99.999% at particle sizes between 0.1 to 0.3µm.	ASTS-030U	ASVLFP536-030U	ASTS-030U (2)
Spill Tray (Polypropylene)	Removable for easy cleaning.	TRAY-P5-24	TRAY-P5-36	TRAY-P5-48
Spill Tray (Stainless Steel)	Removable for easy cleaning.	TRAY-P5-24-SS	TRAY-P5-36-SS	TRAY-P5-48-SS
Base Stand, Mobile, with Casters	Provides a lower storage shelf; accommodates wheelchair access. Locking casters fix the hood in place.	CART-25	CART-36	CART-50
Base Cabinet, Fixed (Metal)	Provides storage space below.	CART-MCC-25	CART-MCC-36	CART-MCC-50
Base Cabinet, Fixed (Polypropylene)	Provides storage space below.	CART-SSC-25	CART-SSC-36	CART-SSC-50
Fire Safety Cabinet Base	Flame resistant safe storage for combustible and flammable liquids	CART-FSC-25	CART-FSC-36	CART-FSC-50
Polypropylene Construction*	Cabinets are available in all polypropylene construction. Contact Air Science for information.	FLOW-24-PP	FLOW-36-PP	FLOW-48-PP

Factory installed; specify when ordering.

Performance & Selection (p.4)

C :(: '. ' / ')

Options and Accessories (p.8

FLOW Laminar Flow Cabinets 24 • 36 • 48 OPTIONS & ACCESSORIES

9

WARRANTY

This product is protected by the Air Science **Legacy Lifetime Warranty™** which starts on the date of shipment from our factory. This limited warranty is the result of thousands of successful Air Science production applications in pharmaceutical, laboratory, forensic, industrial and educational applications.

This warranty covers defects in materials and workmanship. Our liability under the Legacy Lifetime Warranty is, at our option, to repair or replace any defective parts of this equipment if you document that these parts were defective at the time we sold the product to you. Normal conditions apply.



For details visit the <u>Service section</u> of our website at www.airscience.com.

STANDARDS & COMPLIANCE			
Quality Management Systems	ISO 9001:2015		
Environment	ISO 14001:2015 ENERGY STAR® Partner		
Cabinet Performance	IEST-RP-CC002.3 AS 1386.5		
Air Quality	ISO 14644-1, Class 5		
Filtration	IEST-RP-CC034.1 IEST-RP-CC034.3 IEST-RP-CC034.7 EN 1822		
Electrical Safety	UL-C-61010-1 CE Mark RoHS Exempt under EEE Catergory 9		



120 6th Street \ Fort Myers, FL 33907

T. 239-489-0024 \ Toll Free. 800-306-0656 \ F. 800-306-0677

www.airscience.com

The information contained in this manual and the accompanying product are copyrighted and all rights are reserved by Air Science. Air Science reserves the right to make periodic minor design changes without obligation to notify any person or entity of such change.





