

Purair®

FLOW
PCR
LF Series

Consumables Reorder Form

Maintaining your Air Science Laminar Flow Cabinet just became easier.

Four easy ways to reorder service items for your flow cabinet.



1. Fill out the form on the next page and fax it to 800-306-0677 or 239-489-0922; OR
2. Send an email to Air Science: reorder@air-science.com with the part numbers you need to order or attach a purchase order; OR
3. Select the "Please review my history" option on the form and we will review your units history and recommend service items; OR
4. Complete your entire order online at www.filterfilters.com.





CONSUMABLES SELECTION

Option	Part Number	Description	Price	Quantity	Total
--------	-------------	-------------	-------	----------	-------

Purair Flow and PCR Cabinets

	ASTS-FLOW-PRF	Replacement Pre-filters - 6 pack. (Change every 2-3 months)	126.00	x _____	= _____
	ASTS-030	Replacement Hepa Filter. (Change every 12 months)	220.00	x _____	= _____

Purair - LF Series Cabinets

	Prefilters	Replacement Pre-filters - 6 pack. (Change every 2-3 months)			
	VLF72-PRF	Purair - LF Series, VLF-72	126.00	x _____	= _____
	VLF48-PRF	Purair - LF Series, VLF-48	126.00	x _____	= _____
	VLF36-PRF	Purair - LF Series, VLF-36	126.00	x _____	= _____
	HLF96-PRF	Purair - LF Series, HLF-96 (Standard and XT Units)	126.00	x _____	= _____
	HLF72-PRF	Purair - LF Series, HLF-72 (Standard and XT Units)	126.00	x _____	= _____
	HLF48-PRF	Purair - LF Series, HLF-48 (Standard and XT Units)	126.00	x _____	= _____
	HLF36-PRF	Purair - LF Series, HLF-36 (Standard and XT Units)	126.00	x _____	= _____
	Hepa / ULPA Main Filters	Replacement Hepa or ULPA Main Filters (Change every 12 months)			
	ASLF-72-030	Purair - LF Series, VLF-72. ULPA Filter.	695.00	x _____	= _____
	ASLF-48-030	Purair - LF Series, VLF-48. ULPA Filter.	595.00	x _____	= _____
	ASLF-36-030	Purair - LF Series, VLF-36. ULPA Filter.	545.00	x _____	= _____
	ASLF-72-030	Purair - LF Series, HLF-72. ULPA Filter.	695.00	x _____	= _____
	ASLF-48-030	Purair - LF Series, HLF-48 and HLF-96 (2 needed). ULPA Filter.	595.00	x _____	= _____
	ASLF-36-030	Purair - LF Series, HLF-36 ULPA Filter.	545.00	x _____	= _____
	ASLF-72-030-XT	Purair - LF Series, HLF-72XT. ULPA Filter.	750.00	x _____	= _____
	ASLF-48-030-XT	Purair - LF Series, HLF-48XT and HLF-96XT (2 needed). ULPA Filter.	650.00	x _____	= _____
	ASLF-36-030-XT	Purair - LF Series, HLF-36XT. ULPA Filter.	600.00	x _____	= _____

Accessories for all Units

	UVB30	UV Bulbs - 4 pack.	260.00	x _____	= _____
	IV-(model number)	IV-Bar with 'S' Hooks.	149.00	x _____	= _____
	SHOOK	Extra Bag of 'S' Hooks.	36.00	x _____	= _____
	SF	Service Fixtures.	88.00	x _____	= _____

TOTAL: _____

PURCHASE INFORMATION

Purchaser's Name: _____

Phone: _____ Email: _____

Ship to Name: _____

Ship to Address: _____

Credit Card Number (Not needed with PO): _____

Expiration Date: _____

CV Code: _____ Billing Zip Code: _____

Purchase Order Number: _____

Shipping will be charged unless your shipping account number is provided here: _____

If you do not know what you need to order to service your unit and would like us to review your units history, please provide your contact information, company name, and (if easily obtained) the serial number located on the top right-hand side of the unit by the power cord: _____



PO Box 62296 • Fort Myers, FL 33907
T/239.489.0024 • Toll Free/800.306.0656 • F/800.306.0677
www.airscience.com



Schedule
Contact GS077802



Made in the U.S.A.