

PURAIR® BASIC

	List Price	Speed Control	Filter Included	Set Up	Difficulty	Stock	Manpower	Tools Required	Warranty
PURAIR 5-24 XT	£1,254.00	Yes	Yes	15 MIN	EASY	In Stock	2	N/A	2 Years Parts/Labour
P5-24 TRAY	£83.00								
CART P5-24	£350.00								
PURAIR 5-36 XT	£1,464.00	Yes	Yes	15 MIN	EASY	In Stock	2	N/A	2 Years Parts/Labour
P5-36 TRAY	£93.00								
CART P15	£410.00								
PURAIR 5-48 XT	£1,777.00	Yes	Yes	15 MIN	EASY	In Stock	2	N/A	2 Years Parts/Labour
P5-48 TRAY	£110.00								
CART P20	£487.00								

FORMULA

Each Purair® Basic includes one set of our Multiplex impregnated carbon main filters to absorb specific contaminants as indicated. Please select the most appropriate filter formula for your application. For filter selection help or questions, please complete our chemical compatibility form [by clicking here](#).

GP Plus!	Solvents, hydrocarbon vapors	ASTS-001
ACI Plus!	Acids, cyanide, sulphur vapors	ASTS-015
AMM	Ammonia vapors	ASTS-012
FOR	Aldehyde vapors	ASTS-013
HEPA	Hepa	ASTS-030
Others Available		



Purair® 5

PURAIR® ADVANCED

	List Price	Speed Control	Filter Included	Set Up	Difficulty	Stock	Manpower	Tools Required	Warranty
PURAIR 10 W/FSA	£1,595.00	Yes	No	30 MIN	EASY	In Stock	2	✗	2 Years Parts/Labour
CART P10	£375.00								
PURAIR 15 W/FSA	£2,099.00	Yes	No	30 MIN	EASY	In Stock	2	✗	2 Years Parts/Labour
CART P15	£410.00								

FORMULA

It is required that one set of our Multiplex impregnated carbon main filters must be selected to absorb specific contaminants as indicated. Please select the most appropriate filter formula for your application. For filter selection help or questions, please complete our chemical compatibility form [by clicking here](#).

GP Plus!	Solvents, hydrocarbon vapors	ASTM-001	£174.00
ACI Plus!	Acids, cyanide, sulphur vapors	ASTM-015	£236.00
AMM	Ammonia vapors	ASTM-012	£197.00
FOR	Aldehyde vapors	ASTM-013	£247.00
EDU	Educational	ASTM-100	£190.00
HEPA	Hepa	ASTM-030	£260.00
MULTI	GP Plus! and Hepa Multilayer	ASTM-200	£260.00
Others Available			



Purair® 10

21 Days

Many other models as seen in the catalogue are available with similar pricing and typical lead time of 21 days Call or [click here for quote](#).

DWS™

	List Price	Speed Control	Filter Included	Set Up	Difficulty	Stock	Manpower	Tools Required	Warranty
DWS24	£1,777.00	Yes	Yes						2 Years Parts/Labour
CART P5-24	£350.00								

FORMULA

Each DWS™ unit includes one set of our Multiplex impregnated carbon main filters to absorb specific contaminants as indicated. Please select the most appropriate filter formula for your application. For filter selection help or questions, please complete our chemical compatibility form **by clicking here**.

GP Plus!	Solvents, hydrocarbon vapors	AST8-001
ACI Plus!	Acids, cyanide, sulphur vapors	AST8-015
AMM	Ammonia vapors	AST8-012
FOR	Aldehyde vapors	AST8-013
EDU	Educational	AST8-100
HEPA	Hepa	AST8-030
MULTI	Hepa/carbon	AST8-200
Others Available		



DWS24

Many other models as seen in the catalogue are available with similar pricing and typical lead time of 21 days Call or **click here for quote**.

SAFEKEEPER®

	List Price	Speed Control	Filter Included	Set Up	Difficulty	Stock	Manpower	Tools Required	Warranty
FDC-006	£2,573.00	No	Yes						2 Years Parts/Labour
FDC-006XL	£2,995.00	No	Yes						2 Years Parts/Labour

Many other models as seen in the catalogue are available with similar pricing and typical lead time of 21 days Call or **click here for quote**.

PURAIR® FLOW

	List Price	Speed Control	Filter Included	Set Up	Difficulty	Stock	Manpower	Tools Required	Warranty
FLOW-24	£929.00	Yes	Yes					N/A	2 Years Parts/Labour

Many other models as seen in the catalogue are available with similar pricing and typical lead time of 21 days Call or **click here for quote**.



Purair® FLOW-24



Safekeeper® FDC-006XL

PURAIR® LF SERIES

	List Price	Speed Control	Filter Included	Set Up	Difficulty	Stock	Manpower	Tools Required	Warranty
VLF-48	£1,950.00	Yes	Yes	15 MIN	EASY	In Stock	3 people	N/A	2 Years Parts/Labour
BASESTAND VLF-BW-48	£310.00								
UV KIT	£399.00								

21 Days Many other models as seen in the catalogue are available with similar pricing and typical lead time of 21 days Call or [click here for quote](#).


EDU™ MOBILE

	List Price	Speed Control	Filter Included	Set Up	Difficulty	Stock	Manpower	Tools Required	Warranty
EDU-M	£1,999.00	Yes	Yes	15 MIN	EASY	21 Days	2 people	N/A	2 Years Parts/Labour

21 Days Many other models as seen in the catalogue are available with similar pricing and typical lead time of 21 days Call or [click here for quote](#).



VLF-48



EDU-M



Air Science® Technologies Ltd.

Office 4 Attn: Jenny Coffin
 Ocean World House • 1 Paynes Road
 Southampton, SO15 3DL, UK • T/(0) 844.740.1875
 UK@airscience.com

Air Science® USA LLC.

PO Box 62296 • Fort Myers, FL 33907
 T/239.489.0024 • Toll Free/800.306.0656 • F/800.306.0677
 www.airscience.com

