

Consumables Reorder Form

Maintaining your Air Science basic fume just became easier.

Three easy ways to reorder service items for your fume hood.

1. Fill out the form on the next page and fax it to 800-306-0677 or 239-489-0922; OR
2. Send an email to Air Science: reorder@air-science.com with the part numbers you need to order or attach a purchase order; OR
3. Complete your entire order online at www.filtcofilters.com.



Purair P5-24-XT

Purair P5-24-XT (RX)





Purair P5-36-XT (NANO)

Purair P5-24-XT SafeSEARCH



CONSUMABLES SELECTION

Option	Part Number	Description	Price	Quantity	Total
	ASTS-PRF	Replacement Pre-filters - 6 pack. (Change every 2-3 months)	139.00	× _____ = _____	
	ASTS-	Replacement Main Filter. (Change every 12 months - select specific filter below) NOTE: See guide below to determine number of filters required.			
	ASTS-001	GP Filter	248.00	× _____ = _____	
	ASTS-030	HEPA Filter	306.00	× _____ = _____	
	ASTS-015	ACI Filter	285.00	× _____ = _____	
	ASTS-013	FOR Filter	310.00	× _____ = _____	
	ASTS-012	AMM Filter	285.00	× _____ = _____	
	ASTS-002	ACR Carbon Filter	325.00	× _____ = _____	
	ASTS-003	ACM Carbon Filter	325.00	× _____ = _____	
	ASTS-016	SUL Carbon Filter	285.00	× _____ = _____	
			TOTAL:	_____	

HELP GUIDE



M. Main Filter Area: You need one main filter and one pre-filter for each blue filter door. This example shows that one main filter and one pre-filter is needed.

W. Filter Part Number: Look in the filter window to determine the main filter part number. Occasionally, filters are "double stacked," meaning two filters are seen in each blue filter door window.

Purair P5-24-XT

ACCESSORIES FOR ALL UNITS

Gastec Pump™	495.00	× _____ = _____
Direct Reading Gas Detection Tubes - 10 pack. Specify target gas.	90.00	× _____ = _____
Dwyer Airflow Meter	59.00	× _____ = _____

PURCHASE INFORMATION

Purchaser's Name: _____

Phone: _____ Email: _____

Ship to Name: _____

Ship to Address: _____

Credit Card Number (Not needed with PO): _____

Expiration Date: _____

CV Code: _____ Billing Zip Code: _____

Purchase Order Number: _____

Shipping will be charged unless your shipping account number is provided here: _____

If you do not know what you need to order to service your unit and would like us to review your unit's history, please provide your contact information, company name, and (if easily obtained) the serial number located on the top right-hand side of the unit by the power cord: _____



120 6th Street \ Fort Myers, FL 33907
T. 239-489-0024 \ Toll Free. 800-306-0656 \ F. 800-306-0677
www.airscience.com

