

Consumables Reorder Form

Maintaining your Air Science laminar flow cabinet just became easier.

Three easy ways to reorder service items for your laminar flow cabinet.



1. Fill out the form on the next page and fax it to 800-306-0677 or 239-489-0922; OR
2. Send an email to Air Science: reorder@air-science.com with the part numbers you need to order or attach a purchase order; OR
3. Complete your entire order online at www.filtcofilters.com.



CONSUMABLES SELECTION

Option	Part Number	Description	Price	Quantity	Total
--------	-------------	-------------	-------	----------	-------

Purair FLOW and PCR Cabinets

	ASVLF-PRF	Replacement Pre-filters - 6 pack. (Change every 2-3 months).	139.00	× _____	= _____
	HEPA Main Filters	Replacement HEPA Filter. (Change every 24 months).			
	ASTS-030	Purair – FLOW Series, FLOW-24 and FLOW-48 (2 needed). Purair – PCR Series, PCR-24 and PCR-48 (2 needed).	306.00	× _____	= _____
	ASVLP536-030	Purair – FLOW Series, FLOW-36. Purair – PCR Series, PCR-36.	396.00	× _____	= _____

Purair LF Series Cabinets

	Pre-filters	Replacement Pre-filters - 6 pack. (Change every 2-3 months).			
	VLF72-PRF	Purair – LF Series, VLF-72.	139.00	× _____	= _____
	VLF60-PRF	Purair – LF Series, VLF-60.	139.00	× _____	= _____
	VLF48-PRF	Purair – LF Series, VLF-48 and VLF-96.	139.00	× _____	= _____
	VLF36-PRF	Purair – LF Series, VLF-36.	139.00	× _____	= _____
	VLF24-PRF	Purair – LF Series, VLF-24.	139.00	× _____	= _____
	HLF72-PRF	Purair – LF Series, HLF-72 (Standard and XT Units).	139.00	× _____	= _____
	HLF60-PRF	Purair – LF Series, HLF-60 (Standard and XT Units).	139.00	× _____	= _____
	HLF48-PRF	Purair – LF Series, HLF-48 (Standard and XT Units) and HLF-96 (Standard and XT Units) (2 needed).	139.00	× _____	= _____
	HLF36-PRF	Purair – LF Series, HLF-36 (Standard and XT Units).	139.00	× _____	= _____
	HLF24-PRF	Purair – LF Series, HLF-24 (Standard and XT Units).	139.00	× _____	= _____
	HEPA / ULPA Main Filters	Replacement ULPA Main Filters (Change every 24 months).			
	ASLF72-030	Purair – LF Series, VLF-72 and HLF-72. ULPA Filter.	795.00	× _____	= _____
	ASLF60-030	Purair – LF Series, VLF-60 and HLF-60. ULPA Filter.	745.00	× _____	= _____
	ASLF48-030	Purair – LF Series, VLF-48 and VLF-96 (2 needed). ULPA Filter.	695.00	× _____	= _____
	ASLF36-030	Purair – LF Series, VLF-36. ULPA Filter.	645.00	× _____	= _____
	ASLF24-030	Purair – LF Series, VLF-24 and HLF-24. ULPA Filter.	495.00	× _____	= _____
	ASLF72-030-XT	Purair – LF Series, HLF-72XT. ULPA Filter.	850.00	× _____	= _____
	ASLF60-030XT	Purair – LF Series, HLF-60XT. ULPA Filter.	800.00	× _____	= _____
	ASLF48-030-XT	Purair – LF Series, HLF-48XT and HLF-96XT (2 needed). ULPA Filter.	750.00	× _____	= _____
	ASLF36-030-XT	Purair – LF Series, HLF-36XT. ULPA Filter.	700.00	× _____	= _____
	ASLF24-030XT	Purair – LF Series, HLF-24XT. ULPA Filter.	550.00	× _____	= _____

ACCESSORIES FOR ALL UNITS

UVBULB-UV-15 UVBULB-UV-30	UV Bulbs - 4 pack. Provide model #.	260.00 × _____ = _____
IV-(model number)	IV-Bar with 'S' Hooks.	149.00 × _____ = _____
SHOOK	Extra Bag of 'S' Hooks.	36.00 × _____ = _____
SF	Service Fixtures.	88.00 × _____ = _____
		TOTAL: _____

PURCHASE INFORMATION

Purchaser's Name: _____ Phone: _____ Email: _____

Ship to Name: _____ Ship to Address: _____

Credit Card Number (Not needed with PO): _____

Expiration Date: _____ CV Code: _____ Billing Zip Code: _____

Purchase Order Number: _____

Shipping will be charged unless your shipping account number is provided here: _____

If you do not know what you need to order to service your unit and would like us to review your unit's history, please provide your contact information, company name, and (if easily obtained) the serial number located on the top right-hand side of the unit by the power cord: _____



120 6th Street \ Fort Myers, FL 33907
 T. 239-489-0024 \ Toll Free. 800-306-0656 \ F. 800-306-0677
 www.airscience.com

