

# DWS™

## Consumables Reorder Form

Maintaining your Air Science downflow workstation just became easier.

### Three easy ways to reorder service items for your fume hood.


1. Fill out the form on the next page and fax it to 800-306-0677 or 239-489-0922; OR
2. Send an email to Air Science: [reorder@air-science.com](mailto:reorder@air-science.com) with the part numbers you need to order or attach a purchase order; OR
3. Complete your entire order online at [www.filtcofilters.com](http://www.filtcofilters.com).



DWS24



CONSUMABLES SELECTION

Option	Part Number	Description	Price	Quantity	Total
	AST8-PRF	Replacement DWS24 (24" wide unit) Pre-filters - 6 pack. (Change every 2-3 months)	126.00	x _____	= _____
	AST8-PRFXL	Replacement DWS36 (36" wide unit) Pre-filters - 6 pack. (Change every 2-3 months)	174.00	x _____	= _____
	AST8-PRFXXL	Replacement DWS48 (48" wide unit) Pre-filters - 6 pack. (Change every 2-3 months)	174.00	x _____	= _____
	AST8-OR AST8- -XL	Replacement Main Filter. (Change every 12 months - select specific filter below) NOTE: See guide below to determine number of filters required.			
	AST8-001	DWS24 GP Filter	355.00	x _____	= _____
	AST8-001XL	DWS36/48 GP Filter	425.00	x _____	= _____
	AST8-015	DWS24 ACI Filter	425.00	x _____	= _____
	AST8-015XL	DWS36/48 ACI Filter	495.00	x _____	= _____
	AST8-013	DWS24 FOR Filter	445.00	x _____	= _____
	AST8-013XL	DWS36/48 FOR Filter	495.00	x _____	= _____
	AST8-012	DWS24 AMM Filter	425.00	x _____	= _____
	AST8-012XL	DWS36/48 AMM Filter	495.00	x _____	= _____
	AST8-200	DWS24 GP/HEPA Filter	455.00	x _____	= _____
	AST8-200XL	DWS36/48 GP/HEPA Filter	525.00	x _____	= _____
	AST8-030	DWS24 HEPA Filter	445.00	x _____	= _____
	AST8-030XL	DWS36/48 HEPA Filter	495.00	x _____	= _____
	AST8-110	DWS24 Multi-Mix Filter	495.00	x _____	= _____
	AST8-110XL	DWS36/48 Multi-Mix Filter	525.00	x _____	= _____
			TOTAL:	_____	_____

HELP GUIDE



DWS24

**W. Main Filter Area:** DWS units take one main filter and one pre-filter. Look in the filter window to determine the main filter part number.

*NOTE: You will notice that units that are 36 or 48 inches wide (DWS36 or DWS48) have an "XL" at end of the filter part number. 24 inch wide units (DWS24) do not have an "XL" as part of the filter part number.*

**ACCESSORIES**

Gastec Pump™	495.00	x	_____	=	_____
Direct Reading Gas Detection Tubes - 10 pack. Specify target gas.	90.00	x	_____	=	_____
Electronic Filter Saturation Alarm (Boxed unit, sits on top of unit near dirty airstream. 110 - 130v 60Hz or 230v 50Hz)	445.00	x	_____	=	_____
SIDE-DWS - Side Panels to Enhance Containment, set of 2.	175.00	x	_____	=	_____

**PURCHASE INFORMATION**

Purchaser's Name: \_\_\_\_\_ Phone: \_\_\_\_\_ Email: \_\_\_\_\_

Ship to Name: \_\_\_\_\_ Ship to Address: \_\_\_\_\_

Credit Card Number (Not needed with PO): \_\_\_\_\_

Expiration Date: \_\_\_\_\_ CV Code: \_\_\_\_\_ Billing Zip Code: \_\_\_\_\_

Purchase Order Number: \_\_\_\_\_ Shipping will be charged unless your shipping account number is provided here: \_\_\_\_\_

If you do not know what you need to order to service your unit and would like us to review your unit's history, please provide your contact information, company name, and (if easily obtained) the serial number located on the top right-hand side of the unit by the power cord: \_\_\_\_\_



120 6th Street \ Fort Myers, FL 33907  
**T. 239-489-0024 \ Toll Free. 800-306-0656 \ F. 800-306-0677**  
[www.airscience.com](http://www.airscience.com)

