

Drugkeeper™ SafeSTORE™ SafeSTORE™ Mini

Consumables Reorder Form

Maintaining your Air Science special purpose storage cabinet just became easier.

Three easy ways to reorder service items for your fume hood.

1. Fill out the form on the next page and fax it to 800-306-0677 or 239-489-0922; OR
2. Send an email to Air Science: reorder@air-science.com with the part numbers you need to order or attach a purchase order; OR
3. Complete your entire order online at www.filtcofilters.com.







Drugkeeper 64T

SafeSTORE-34T



SafeSTORE Mini



CONSUMABLES SELECTION

Option	Part Number	Description	Price	Quantity	Total
	SS-PRF	Replacement Door Pre-filters for 34T - 6 pack. (Change every 2-3 months)	139.00	× _____	= _____
	SS64T-PRF	Replacement Door Pre-filters for 64T - 6 pack. (Change every 2-3 months)	139.00	× _____	= _____
	ASTS-PRF	Replacement Interior Pre-filters - 6 pack. (Change every 2-3 months)	139.00	× _____	= _____
	ASTS-	Replacement Safety Filter.(Chage every 12 months - select specific filter below) NOTE: See guide on following page to determine number of filters required.			
	ASTS-001	GP Filter	248.00	× _____	= _____
	ASTS-012	AMM Filter	285.00	× _____	= _____
	ASTS-013	FOR Filter	310.00	× _____	= _____
	ASTS-015	ACI Filter	285.00	× _____	= _____
	ASTS-030	HEPA Filter	306.00	× _____	= _____

The following filters are for the MINI units only:

	AP60-PRF	Replacement Interior Pre-filters - 6 pack. (Change every 2-3 months)	98.00	× _____	= _____
	AP60-	Replacement Safety Filter.(Chage every 12 months - select specific filter below) NOTE: See guide on following page to determine number of filters required.			
	AP60-001	GP Filter	145.00	× _____	= _____
	AP60-013	FOR Filter	168.00	× _____	= _____
	AP60-015	ACI Filter	155.00	× _____	= _____
	AP60-030	HEPA Filter	264.00	× _____	= _____

TOTAL: _____

HELP GUIDE



M. Main Filter Area: You need one HEPA filter and one carbon filter for each blue filter door. This unit happens to have two blue filter doors. You would need two carbon filters, two HEPA filters, and two pre-filters for this main filter area.

D. Door Filter Area: One pre-filter is needed for each door filter. This unit happens to need two pre-filters.

Units produced before 2005 take a different filter (part number ASTM-200). If you think your unit may take this filter, let us know. We will check your units history so that the correct filter is shipped.

Drugkeeper 64T

DOES MY UNIT HAVE THE OPTIONAL SAFETY FILTER?

Look at the top of the unit. If you can see the filter bolted by screws, your unit has this option.



In the photo above, the filter can be readily seen on the top of the unit. If you just see solid white on top, it does not have the optional safety filter.

PURCHASE INFORMATION

Purchaser's Name: _____ Phone: _____ Email: _____

Ship to Name: _____ Ship to Address: _____

Credit Card Number (Not needed with PO): _____

Expiration Date: _____ CV Code: _____ Billing Zip Code: _____

Purchase Order Number: _____ Shipping will be charged unless your shipping account number is provided here: _____

If you do not know what you need to order to service your unit and would like us to review your unit's history, please provide your contact information, company name, and (if easily obtained) the serial number located on the top right-hand side of the unit by the power cord: _____



120 6th Street \ Fort Myers, FL 33907
 T. 239-489-0024 \ Toll Free. 800-306-0656 \ F. 800-306-0677
 www.airscience.com

