

## Consumables Reorder Form

Maintaining your Air Science Laminar Flow Cabinet just became easier.

Four easy ways to reorder service items for your laminar flow cabinet.

1. Fill out the form on the next page and fax it to 800-306-0677 or 239-489-0922; OR
2. Send an email to Air Science: [reorder@air-science.com](mailto:reorder@air-science.com) with the part numbers you need to order or attach a purchase order; OR
3. Select the "Please review my history" option on the form and we will review your units history and recommend service items; OR
4. Complete your entire order online at [www.filtcofilters.com](http://www.filtcofilters.com).





Purair® FLOW-24

Purair PCR-24



VLF-48



CONSUMABLES SELECTION

Option	Part Number	Description	Price	Quantity	Total
<b>Purair FLOW and PCR Cabinets</b>					
	ASTS-FLOW-PRF	Replacement Pre-filters - 6 pack. (Change every 2-3 months).	126.00	× _____	= _____
	Hepa / ULPA Main Filters	Replacement Hepa Filter. (Change every 12 months).			
	ASTS-030	Purair – FLOW Series, FLOW-24 and FLOW-48 (2 needed). Purair – PCR Series, PCR-24 and PCR-48 (2 needed).	220.00	× _____	= _____
	ASVLP536-030	Purair – FLOW Series, FLOW-36. Purair – PCR Series, PCR-36.	298.00	× _____	= _____

**Purair LF Series Cabinets**

	Prefilters	Replacement Pre-filters - 6 pack. (Change every 2-3 months).			
	ASVLF-PRF	Purair – LF Series, VLF.	126.00	× _____	= _____
	VLF72-PRF	Purair – LF Series, VLF-72.	126.00	× _____	= _____
	VLF60-PRF	Purair – LF Series, VLF-60.	126.00	× _____	= _____
	VLF48-PRF	Purair – LF Series, VLF-48 and VLF-96.	126.00	× _____	= _____
	VLF36-PRF	Purair – LF Series, VLF-36.	126.00	× _____	= _____
	VLF24-PRF	Purair – LF Series, VLF-24.	126.00	× _____	= _____
	HLF96-PRF	Purair – LF Series, HLF-96 (Standard and XT Units).	126.00	× _____	= _____
	HLF72-PRF	Purair – LF Series, HLF-72 (Standard and XT Units).	126.00	× _____	= _____
	HLF60-PRF	Purair – LF Series, HLF-60 (Standard and XT Units).	126.00	× _____	= _____
	HLF48-PRF	Purair – LF Series, HLF-48 (Standard and XT Units).	126.00	× _____	= _____
	HLF36-PRF	Purair – LF Series, HLF-36 (Standard and XT Units).	126.00	× _____	= _____
	HLF24-PRF	Purair – LF Series, HLF-24 (Standard and XT Units).	126.00	× _____	= _____
	Hepa / ULPA Main Filters	Replacement Hepa or ULPA Main Filters (Change every 12 months).			
	ASLF-72-030	Purair – LF Series, VLF-72 and HLF-72. ULPA Filter.	695.00	× _____	= _____
	ASLF-60-030	Purair – LF Series, VLF-60 and HLF-60. ULPA Filter.	650.00	× _____	= _____
	ASLF-48-030	Purair – LF Series, VLF-48 and VLF-96 (2 needed). ULPA Filter.	595.00	× _____	= _____
	ASLF-36-030	Purair – LF Series, VLF-36. ULPA Filter.	545.00	× _____	= _____
	ASLF-24-030	Purair – LF Series, VLF-24 and HLF-24. ULPA Filter.	395.00	× _____	= _____
	ASLF-48-030	Purair – LF Series, HLF-48 and HLF-96 (2 needed). ULPA Filter.	595.00	× _____	= _____
	ASLF-36-030	Purair – LF Series, HLF-36 ULPA Filter.	545.00	× _____	= _____
	ASLF-72-030-XT	Purair – LF Series, HLF-72XT. ULPA Filter.	750.00	× _____	= _____
	ASLF-60-030XT	Purair – LF Series, HLF-60XT. ULPA Filter.	695.00	× _____	= _____
	ASLF-48-030-XT	Purair – LF Series, HLF-48XT and HLF-96XT (2 needed). ULPA Filter.	650.00	× _____	= _____
	ASLF-36-030-XT	Purair – LF Series, HLF-36XT. ULPA Filter.	600.00	× _____	= _____
	ASLF-24-030XT	Purair – LF Series, HLF-24XT. ULPA Filter.	395.00	× _____	= _____

ACCESSORIES FOR ALL UNITS

UVB30	UV Bulbs - 4 pack.	260.00	×	_____	=	_____
IV-(model number)	IV-Bar with 'S' Hooks.	149.00	×	_____	=	_____
SHOOK	Extra Bag of 'S' Hooks.	36.00	×	_____	=	_____
SF	Service Fixtures.	88.00	×	_____	=	_____
<b>TOTAL:</b>						_____

PURCHASE INFORMATION

Purchaser's Name: \_\_\_\_\_ Phone: \_\_\_\_\_ Email: \_\_\_\_\_

Ship to Name: \_\_\_\_\_ Ship to Address: \_\_\_\_\_

Credit Card Number (Not needed with PO): \_\_\_\_\_

Expiration Date: \_\_\_\_\_ CV Code: \_\_\_\_\_ Billing Zip Code: \_\_\_\_\_

Purchase Order Number: \_\_\_\_\_

Shipping will be charged unless your shipping account number is provided here: \_\_\_\_\_

If you do not know what you need to order to service your unit and would like us to review your unit's history, please provide your contact information, company name, and (if easily obtained) the serial number located on the top right-hand side of the unit by the power cord: \_\_\_\_\_



120 6th Street \ Fort Myers, FL 33907  
 T. 239-489-0024 \ Toll Free. 800-306-0656 \ F. 800-306-0677  
 www.airscience.com

