

# SafeFume™

## Consumables Reorder Form

Maintaining your Air Science Automatic Cyanoacrylate Fuming Chamber just became easier.

### Four easy ways to reorder service items for your fuming chamber.

1. Fill out the form on the next page and fax it to 800-306-0677 or 239-489-0922; OR
2. Send an email to Air Science: [reorder@air-science.com](mailto:reorder@air-science.com) with the part numbers you need to order or attach a purchase order; OR
3. Select the "Please review my history" option on the form and we will review your units history and recommend service items; OR
4. Complete your entire order online at [www.filtcofilters.com](http://www.filtcofilters.com).














SafeFume™ CA30S

SafeFume™ 360  
ARV-24S



**CONSUMABLES SELECTION**

Option	Part Number	Description	Price	Quantity	Total
	ASTM-001	Replacement Carbon Filter standard Safefume (CA) models excluding CA24S and Safefume 360 (ARV) models. (Change every 12 months)	355.00	× _____ = _____	
	ASCAMINI-001	Replacement Carbon Filter for CA24S. (Change every 12 months)	138.00	× _____ = _____	
	AP60H-001	Replacement Carbon Filter for Safefume 360 (ARV) models excluding ARV-72XL. (Change every 12 months)	126.00	× _____ = _____	
	ASTS-001	Replacement Carbon Filter for ARV-72XL. (Change every 12 months)	220.00	× _____ = _____	
	ASTM-PRF	Replacement Pre-filters for standard Safefume (CA) models excluding Safefume 360 (ARV) models. - 6 pack. (Change every 2-3 months)	174.00	× _____ = _____	
	AP60-PRF	Replacement Pre-filters for CA24S and Safefume 360 (ARV) models excluding ARV-72XL. - 6 pack. (Change every 2-3 months)	90.00	× _____ = _____	
	CA-PRF	Replacement Circulation Fan Pre-filters - 10 pack. (Change every 2-3 months)	88.00	× _____ = _____	
	72452	57 mm Disposable Aluminum Fuming Dishes - pack of 100.	12.00	× _____ = _____	
	AA-490	Cyanoacrylate - 20 g bottle. Use our super thin glue and you will see quick results.	7.00	× _____ = _____	
	AA-500	Cyanoacrylate - 500 g bottle.	85.00	× _____ = _____	
	AA-514	Cyanoacrylate - empty 50 g bottles for dispensing from larger bottle.	3.00	× _____ = _____	
	WF2	Wicking Filter. (Located in humidifier)	16.00	× _____ = _____	
	8442	Evidence Hanging Clips. (6 per pack)	9.00	× _____ = _____	
	6171844	Versatile clip and dry system uses 26 high-quality firm-grip clips that hang from strong white chain-like extension from rod.	26.00	× _____ = _____	
	WESRX63080D	Pre-treated Wipes - pack of 25. Specially formulated to protect new, unused fuming chambers constructed of polycarbonate material or glass. Wipe the pad across all interior glass or polycarbonate surfaces, allow it to dry and apply a second coat. When fuming is complete, wipe the interior with damp cloth.	25.00	× _____ = _____	
	PP-2OZ	Plastic Cleaning Starter Kit 2 polish cloths + Polish 1 - Clean & Shine - 2 oz / Polish 2 - Fine Scratch Remover - 2 oz Polish 3 - Heavy Scratch Remover - 2 oz / Two polish mates.	36.00	× _____ = _____	
	HOTPLATE-110	Replacement Hot Plate	60.00	× _____ = _____	
	P10XL-CART	CA24S/CA30S Heavy Duty Steel Cart with Lockable Casters	595.00	× _____ = _____	
	P20-CART	CA48S Heavy Duty Steel Cart with Lockable Casters	695.00	× _____ = _____	
	P25-CART	CA60S Heavy Duty Steel Cart with Lockable Casters	735.00	× _____ = _____	
			TOTAL:	_____	

PURCHASE INFORMATION

Purchaser's Name: \_\_\_\_\_

Phone: \_\_\_\_\_ Email: \_\_\_\_\_

Ship to Name: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

Credit Card Number (Not needed with PO): \_\_\_\_\_

\_\_\_\_\_

Expiration Date: \_\_\_\_\_

CV Code: \_\_\_\_\_ Billing Zip Code: \_\_\_\_\_

Purchase Order Number: \_\_\_\_\_

Shipping will be charged unless your shipping account number is provided here: \_\_\_\_\_

If you do not know what you need to order to service your unit and would like us to review your units history, please provide your contact information, company name, and (if easily obtained) the serial number located on the top right-hand side of the unit by the power cord: \_\_\_\_\_



PO Box 62296 • Fort Myers, FL 33907  
T/239.489.0024 • Toll Free/800.306.0656 • F/800.306.0677  
www.airscience.com



Schedule Contract GS-07F5832P

