



**NOTES :**

- 1.0 CUSTOMER
- 2.0 SUBMIT DRAWING FOR APPROVAL.
- 3.0 STANDARD FEATURES INCLUDED.
- 3.1 BODY CONSTRUCTION, WELDED POLYPRO.
- 3.2 POWER SUPPLY 120 VOLT 60 HZ.
- 3.3 ON / OFF SWITCH ON CONTROL PANEL AND SATURATION FILTER ALARM.
- 3.4 STANDARD ACCESSORIES.



MAIN FILTER

A-A ( 1 : 4 )


TYPICAL  
2" Ø PIPE  
CONNECTION

BLIND  
FLANGE

CONTROL  
PANEL

FAN  
INSIDE

POWER  
INLET



B-B ( 1 : 4 )

Designed by CFG	Checked by	Approved by	Date	Date 2/20/2019	
AIRB SCIENCE USA			AP60 VENT BOX		
			AVB60-Unit		Edition Sheet 1 / 1