





NOTES :

1.0 CUSTOMER:

- 2.0 SUBMIT DRAWINGS FOR APPROVAL.
- 3.0 STANDARD FEATURES INCLUDE:
- 3.1 CABINET BUILT WITH ALUMINUM EXTRUDE WATERTIGHT INSIDE.
- 3.2 FRONT, TOP SIDES AND DOOR WITH TEMPERED GLASS.
- 3.3 CARBON MAIN FILTER.
- 3.4 POWER SUPPLY 120 V 60 HZ.
- 3.5 HOT PLATE ACCELERATOR.
- 3.6 NEBULIZER / HUMIDIFIER.
- 3.7 STANDARD ACCESORIES APPLIED.
- 3.8 MICROPROCESSOR SCREEN TOUCH PAD CONTROL PANEL.
- 4.0 W/O



A-A (1 : 14)

Designed by CFG	Checked by	Approved by	Date 6/19/2019	Date
AIR SCIENCE USA		ARV-60T		
ARV60-G.Assy		Edition	Sheet 1 / 1	