


A-A (1 : 10)



- NOTES :
- 1.0 CUSTOMER:
 - 2.0 SUBMIT DRAWINGS FOR APPROVAL.
 - 3.0 BODY FABRICATED IN WELDED POLYPROPYLENE.
 - 4.0 POWER INLET 110 VOLT, 60 HZ.

Designed by CG	Checked by	Approved by	Date 12/8/2014	Date
AIR SCIENCE USA		DFC 36" DOWNFLOW DRYING CABINET		
		DFDC36Assem	Edition R3	Sheet 1 / 1