



REVISION HISTORY

ZONE	REV	DESCRIPTIO	DATE	APPROVED
DESIGN	1	UPDATE	6/29/2017	P CH




NOTES :

- 1.0 CUSTOMER:
- 2.0 SUBMIT DRAWINGS FOR APPROVAL.
- 3.0 INCLUDED:
- 3.1 BODY BUILT WITH 3/8" WHITE POLYPROPYLENE, SEAMLESS, NON-POROUS WATERTIGHT INSIDE, TERMALLY FUSE.
- 3.2 SEAMLESS POLYPRO WORKSURFACE.
- 3.3 INTERNAL LIGHTING.
- 3.4 SPEED CONTROL.
- 3.5 POWER SWITCH ON/OFF.
- 4.0 ELECTRICAL REQUIREMENT 110 VOLT, 60 HZ.
- 5.0 AST8 015 ACID FILTER.

A-A (1 : 7)

B-B (1 : 7)



Designed by CFG	Checked by	Approved by	Date 6/7/2017	Date
AIR SCIENCE USA		DOWN FLOW WORKSTATION		
DWS.PP24-Gassem			Edition 1	Sheet 1 / 1