



REVISION HISTORY				
ZONE	REV	DESCRIPTION	DATE	APPROVED
FRONT PANEL	1	CONTROL PANEL POSITION	10/04/2016	P CH



HANGING RODS




- NOTES :
- 1.0 CUSTOMER:
  - 2.0 SUBMIT DRAWINGS FOR APPROVAL.
  - 3.0 BODY UNIT FABRICATED IN WELDED POLYPROPYLENE.
  - 4.0 POWER INLET 110 VOLT, 60 HZ.

DOOR OPEN 230 X 360(2)

CLAMP ACCESS FRONT DOORS

CONTROL PANEL(2)



MAIN FILTERS(2)

FRONT METAL DOORS WITH TEMPERATE GLASSES (2)

PRE-FILTER DOORS

SHELVES(4)

FRONT HANDLES(2)

DOOR SERVICES

MAIN FILTERS METAL DOOR

POWER INLETS(2)

FRONT LOCKING CASTORS

Designed by CG	Checked by	Approved by	Date 8/4/2016	Date
AIR SCIENCE USA			FDC-60D DOWNFLOW DRYING CABINET	
DFDC60-GAssem			Edition 1	Sheet 1 / 1