



B-B ( 1 : 5 )



A-A ( 1 : 5 )

- NOTES:**
- A.0 STANDARD ON ALL UNITS:
  - 1.0 CONSTRUCTION: EPOXY COATED MILD STEEL HEAD UNIT AND FRAME.
  - 2.0 FRONT, SIDES AND BACK PANELS MADE OF 1/4" CLEAR POLYCARBONATE.
  - 3.0 MAIN SWITCH ON/OFF.
  - 4.0 CENTRIFUGAL FAN.
  - 5.0 POWER SUPPLY:
  - 5.1 FOR 120 VOLT 60 HZ.
  - 5.2 FOR 230 VOLT 50 HZ.
  - B.0 MAIN FILTER IS APPLICATION-SPECIFIC; MUST ORDER SEPARATELY.



Designed by CFG	Checked by	Approved by	Date 6/25/2019	Date
AIR SCIENCE USA		MINI DOUBLE SAFESTORE		
MiniDouble-SSStg02.GAssy			Edition	Sheet 1 / 1