



NOTES:

- A.0 STANDARD ON ALL UNITS:
- 1.0 CONSTRUCTION: EPOXY COATED MILD STEEL HEAD UNIT AND FRAME.
- 2.0 FRONT, SIDES AND BACK PANELS MADE OF 1/4" CLEAR POLYCARBONATE.
- 3.0 MAIN SWITCH ON/OFF.
- 4.0 CENTRIFUGAL FAN.
- 5.0 POWER SUPPLY:
- 5.1 FOR 120 VOLT 60 HZ.
- 5.2 FOR 230 VOLT 50 HZ.
- B.0 MAIN FILTER IS APPLICATION-SPECIFIC; MUST ORDER SEPARATELY.



A-A (1 : 5)



B-B (1 : 5)

Designed by CF	Checked by	Approved by	Date 6/25/2019	Date	
AIR SCIENCE USA		MINI SINGLE SAFESTORE CSC		Edition	Sheet 1 / 1
		Mini-Stg01.Gassy			