



- NOTES:**
- STANDARD ON ALL UNITS:
  - 1.0 CONSTRUCTION: EPOXY COATED MILD STEEL HEAD UNIT AND COLUMNS.
  - 2.0 FRONT, SIDES AND BACK PANELS MADE ON CLEAR ACRYLIC OF 1/4".
  - 3.0 CENTRIFUGAL FAN
  - 4.0 POWER SUPPLY 110V 60Hz
  - 5.0 SWITCHES, MAIN ON/OFF
  - 6.0 MONITORING, LOW AIRFLOW

Designed by CFG	Checked by	Approved by	Date	Date 4/26/2013	
AIR SCIENCE USA		P5-36SL UNIT			
		P5.P36-MainAssem-R1	Edition	Sheet 1 / 1	