




REVISION HISTORY				
ZONE	REV	DESCRIPTION	DATE	APPROVED
SIDE PANELS	1	TEMPERED GLASSES	5/16/2017	P CH




NOTES:
 CUSTOMER:
 SUBMIT DRAWING FOR APPROVAL

STANDARD ON THIS UNIT:
 1.0 CONSTRUCTION: EPOXY COATED MILD STEEL SIDES PANELS AND HEAD UNIT.
 2.0 FRONT, SIDES AND MADE ON CLEAR PANELS OF 1/4", PP WHITE BACK WALL.
 3.0 TEMPERED GLASSES ON SIDE PANELS.
 4.0 CENTRIFUGAL FANS.
 5.0 POWER SUPPLY 110V 60Hz.
 6.0 MAIN SWITCH ON/OFF.
 7.0 MONITORING, LOW AIRFLOW.



B-B (1 : 10)



A-A (1 : 10)

Designed by CFG	Checked by	Approved by	Date 12/31/2008	Date
AIR SCIENCE USA			PURAIR 40UL DUCTLESS FUME HOOD	
P40UL.-Assem			Edition 1	Sheet 1 / 1