



A-A (1 : 8)



- NOTES:**
- STANDARD ON ALL UNITS:
 - 1.0 CONSTRUCTION: EPOXY COATED MILD STEEL HEAD UNIT AND COLUMNS.
 - 2.0 FRONT TOP, SIDES AND BACK PANELS MADE ON CLEAR ACRYLIC OF 1/4".
 - 3.0 CENTRIFUGAL FAN
 - 4.0 POWER SUPPLY 110V 60Hz
 - 5.0 SWITCHES, MAIN ON/OFF

Designed by CFG	Checked by	Approved by	Date 8/7/2012	Date
AIR SCIENCE USA		PURAIR FLOW PF-24		
		PF-24	Edition	Sheet 1 / 1