




REVISION HISTORY				
ZONE	REV	DESCRIPTION	DATE	APPROVED
SIDE PANELS	1	TEMPERED GLASSES	5/12/2017	P CH



B-B (1 : 10)



A-A (1 : 10)



- NOTES:**
- 1.0 CONSTRUCTION: EPOXY COATED MID STEEL SIDES PANELS AND HEAD UNIT.
 - 2.0 FRONT AND BACK PANELS MADE ON CLEAR ACRYLIC OF 1/4".
 - 3.0 TEMPERED GLASSES ON SIDE PANELS.
 - 4.0 CENTRIFUGAL FAN
 - 5.0 POWER SUPPLY 110V 60Hz
 - 6.0 MAIN SWITCH ON/OFF
 - 7.0 MONITORING, LOW AIRFLOW

Designed by CFG	Checked by	Approved by	Date 9/21/2012	Date
AIR SCIENCE USA		PURAIR 10XT.UL DUCTLESS FUME HOOD		
		Purair10XT-UL.Assem	Edition 1	Sheet 1 / 1