




REVISION HISTORY				
ZONE	REV	DESCRIPTION	DATE	APPROVED
SIDE PANELS	1	TEMPERED GLASSES	5/18/2017	P CH



- NOTES:**
- 1.0 CONSTRUCTION: EPOXY COATED MID STEEL SIDES PANELS AND HEAD UNIT.
 - 2.0 FRONT, SIDES AND BACK PANELS MADE ON CLEAR ACRYLIC OF 1/4".
 - 3.0 TEMPERED GLASSES ON SIDE PANELS.
 - 4.0 CENTRIFUGAL FAN.
 - 5.0 POWER SUPPLY 110V 60Hz.
 - 6.0 MAIN SWITCH ON/OFF.
 - 7.0 MONITORING, LOW AIRFLOW.

Designed by CFG	Checked by	Approved by	Date 9/18/2012	Date
AIR SCIENCE USA		PURAIR 20XT.UL DUCTLESS FUME HOOD		
Purair20XT-UL			Edition 1	Sheet 1 / 1