


- NOTES:**
- STANDARD ON ALL UNITS:
  - 1.0 CONSTRUCTION: EPOXY COATED MID STEEL SIDES PANELS AND HEAD UNIT.
  - 2.0 FRONT, SIDES AND BACK PANELS MADE ON CLEAR ACRYLIC OF 1/4".
  - 3.0 CENTRIFUGAL FAN
  - 4.0 POWER SUPPLY 110V 60Hz
  - 5.0 SWITCHES, MAIN ON/OFF
  - 6.0 MONITORING, LOW AIRFLOW

- CENTRIFUGAL FANS (2)
- FILTERS CLAMPS (2) SYSTEMS
- LAMP MOUNT (2)
- SAFETY FILTERS (2) (OPTIONAL)



Designed by CFG	Checked by	Approved by	Date 12/10/2008	Date
AIR SCIENCE USA		PURAIR 25XT DUCTLESS FUME HOOD		
		Purair25-UL.Assem	Edition 1	Sheet 1 / 1