

The Authority in HEPA support

## S & M Electronics Ltd.

105 Eastcote Road, Pinner,  
Middlesex HA5 1ET, UK.

Tel : (020) 8429 1222

Fax: (020) 8248 4313

info@sandmelectronics.com

# AEROSOL GENERATOR CERTIFICATE OF COMPLIANCE

Certificate No: 20756

Test Date: 06<sup>th</sup> January 2026

Manufacturer: A T I

Model: ATI 6D

Serial No: 34879

Test Due: 06<sup>th</sup> January 2027

Customer: Air Science Technologies Ltd  
Cross farm, Maiscough Lane  
Downholland  
LANCASHIRE L39 7HT

Contact: Mr. Michael Crilly  
Paul

Upon receipt by S & M Electronics Limited, on first inspection & testing, the Generator Air Pressure at laskin Nozel is 20 PSI (Tolerance +/- 2.0 %) & the instrument was found to be in good working order.

S & M Electronics Limited has serviced this instrument in accordance with the Manufacturer's requirements, having been fully tested and found to be in good order.

This equipment has been calibrated using instruments traceable to National Standard. The measurements were correct at the time of calibration.

The calibration was performed in the service area at a room temperature of  $21^{\circ}\text{C} \pm 5^{\circ}\text{C}$ .  
Environmental considerations are not relevant to the performance of the instrument.

### Results:


Compressor Air pressure at Laskin Nozzle 20 P.S.I. (+/- 2%).

NOTE: Indicated reading on pressure gauge of 19.5PSI = ACTUAL 20PSI

Only **Shell Ondina X420 Oil** is recommended for this equipment.

Calibration Completed by: .....*Atul Patel*..... Mr. A. Patel (Engineer)

**UKAS Certified Test Equipment:** Budenberg 100 PSI 6" Pressure Gauge s/n: 58197001 Certificate No: U33939



Website: [www.sandmelectronics.com](http://www.sandmelectronics.com)